




İNOVATÜRK®

THE İNOVATÜRK® STERILE CONTAINER SYSTEM

THE İNOVATÜRK® STERILE CONTAINER SYSTEM



The right INSTRUMENTS

At the right TIME

In sterile CONDITION

At the right PLACE OF USE

At a high level of PROFITABILITY

In a reliable PACKAGING SYSTEM




Fast, reliable and
efficiently organized

**THE
İNOVATÜRK®STERILE
CONTAINER
SYSTEM**

- Integrated permanent germ barrier system made of PTFE designed for 5,000 sterilisation cycles. No filter changes required - saves time and processing costs
- Stable, robust aluminium lid - reliable and long-lasting
- Outstanding cleaning results
- Good drying properties
- Available in five colours (red, blue, green, gold, silver) to help you organize your instrument sets





Article number	Size / colour	Dimensions mm (length x width x height)
JP101	1/1 Lid, red	588 x 281 x 36
JP102	1/1 Lid, blue	588 x 281 x 36
JP103	1/1 Lid, green	588 x 281 x 36
JP104	1/1 Lid, gold	588 x 281 x 36
JP105	1/1 Lid, silver	588 x 281 x 36
JP111	3/4 Lid, red	465 x 281 x 36
JP112	3/4 Lid, blue	465 x 281 x 36
JP113	3/4 Lid, green	465 x 281 x 36
JP114	3/4 Lid, gold	465 x 281 x 36
JP115	3/4 Lid, silver	465 x 281 x 36
JP121	1/2 Lid, red	298 x 281 x 36
JP122	1/2 Lid, blue	298 x 281 x 36
JP123	1/2 Lid, green	298 x 281 x 36
JP124	1/2 Lid, gold	298 x 281 x 36
JP125	1/2 Lid, silver	298 x 281 x 36
JP131	1/2 Lid for JK339, red	298 x 281 x 36
JP132	1/2 Lid for JK339, blue	298 x 281 x 36
JP133	1/2 Lid for JK339, green	298 x 281 x 36
JP134	1/2 Lid for JK339, gold	298 x 281 x 36
JP135	1/2 Lid for JK339, silver	298 x 281 x 36

Compatible with all bottom parts (1/1, 1/2, 3/4) manufactured since 1988
 Consumables and accessories can be found in catalogue no. C40402



Fast, reliable and
efficiently organized
THE INOVATÜRK®
STERILE CONTAINER
SYSTEM

inovatürk



+90 (212) 210 20 20



info@inovaturk.com.tr



www.inovaturk.com.tr

Ataşehir bulvarı, Atatürk Mah. Evinpark Rezidans No:23 a-1 blok d.24 Ataşehir / İSTANBUL