


# LIFECUBE®

RENEW YOUR BODY



# BROCHURE



Content:

1. Cryotherapy.....	3
2. Application area and methods.....	4
3. Advantages of LifeCube.....	5
4. Different LifeCube models.....	6
5. Your personal LifeCube.....	7
6. Advantages of Cosmetic-Solution.....	8



+90 (212) 210 20 20



info@inovaturk.com.tr



www.inovaturk.com.tr


Ataşehir bulvarı, Atatürk Mah. Evinpark Rezidans No:23 a-1 blok d.24 Ataşehir / İSTANBUL

❄️ *Renew your body with the power of freeze* ❄️

Cryotherapy is becoming more popular and successful every day. Even celebrities and athletes are counting on cryotherapy nowadays. This famous method effects body and mind and causes many positive changes. And the best part is that everyone can benefit from it. Due to versatile application possibilities, cryo chambers can be used for medical, fitness, wellness, and beauty purposes. In this context, they can be found in a variety of different places, such as rheumatism clinics and health clinics, wellness centers and spas, sport centers and intensive training centers, beauty salons and cosmetic studios, physiotherapies, and of course with private customers.

normal body temperature

body temperature inside the LifeCube at -80° C



❄️ *There are no reasons to fear the cold. But many to love it.* ❄️

Extreme cold for only a few minutes is neither painful nor unpleasant, but rather enjoyable. The key for a successful cryo therapy is here time management: As long as the customer is only about three minutes in contact with the low temperatures, they are not harmful, but positively effecting their body. Furthermore, as the air inside the chamber is very, very dry, the temperatures do not seem to be that freezing cold. In fact, many customers argue that dives into a swimming pool after a sauna session feel colder than a full body cryo therapy at -110°C.

# APPLICATION AREA AND METHODS

You want to profit from all the benefits cryotherapy has to offer? A short visit is enough to achieve first successes - only 3 minutes inside the cryo chamber are sufficient. Just put on the provided protective cloth, enter the LifeCube and enjoy the pleasant cold.

The use of our LifeCube is easy for the customer as well as the employee. To guarantee operator convenience, we programmed an interactive, almost self-explanatory app. It can be used to regulate the temperature, start the treatment or check the relevant data. Thus, you can simply control the entire cold therapy via app on our supplied tablet or smartwatch.



## SPORT

- sport studio
- fitness studio
- intensive training center
- sports club
- sport organization
- cryo center



## HEALTH


- pain center
- rehab facility
- clinic
- physiotherapy
- ergotherapy
- medical practice



## WELLNESS

- cosmetic studio
- spa hotel
- therme
- swimming pool
- solarium
- sauna

# ADVANTAGES OF LIFECUBE



*A nitrogen cryo chamber is not comparable to an electrical cryo chamber*

Not every whole body cryo chamber is as efficient and effective as our LifeCube. It has the decisive advantage of being an electrical cryo chamber and therefore not being dependent on external media to fully function. Thus, in comparison to nitrogen chambers, our LifeCube succeeds in every way.



## EXTREME TEMPERATURES

LifeCubes can be cooled down to -110°C!



## SAFE APPLICATION

Electrical cryo chambers pose no health risks such as suffocations or frostbites as do nitrogen chambers



## LOW OPERATIONAL COSTS

LifeCubes help you save costs in the long turn with a low power consumption and small maintenance costs



## MADE IN GERMANY

Our cryo chambers are high-class products with the best quality directly from Germany



## INNOVATIVE APP

LifeCubes are practical and easy to operate via interactive app with the supplied tablet and smartwatch

There is not only one LifeCube – there are four of them! Depending on the space you can provide, and the temperatures you want to provide, you can choose between a one-chamber-system, a one-chamber-system with pre-chamber, a two-chamber-system, and a three-chamber-system. While smaller LifeCubes have the advantage of fitting into every room, bigger LifeCubes have the advantage of possessing multiple chambers, which can be tempered differently. This way customers can slowly acclimatize themselves and move from -10°C to -110°C step by step.



## 1 Chamber System

- Temperatures up to -80°C
- Measurements (W x H x D): 100 x 230 x 100 cm
- Weight: 500 kg
- Power consumption: 6 kW

➔ Also available with a pre-chamber!



## 2 Chamber System

- Temperatures up to -110°C
- Measurements (W x H x D): 220 x 230 x 110 cm
- Weight: 700 kg
- Power consumption: 10,5 kW



## 3 Chamber System


- Temperatures up to -110°C
- Measurements (W x H x D): 330 x 230 x 110 cm
- Weight: 1000 kg
- Power consumption: 10,5 kW

# YOUR PERSONAL LIFECUBE

**inovatürk**


In addition to our four standard models, we also offer custom-made LifeCubes. We would love to fulfill your individual wishes and requirements and create your own individual LifeCube!

Our technical drawers will gladly create a draft for you, on which you can see how the LifeCube can be integrated into your studio. Just send us the measurements of your space, explain to us where the LifeCube should stay, and we will give you a glimpse at your future. As the measurements of our LifeCube can be flexible adjusted, corners, less space or a pitched roof are no longer a problem – our electrical cryo chamber fits in every room.



In addition, you can choose between our optional equipment and complement your LifeCube with an intercom system, a bigger aggregate, a video surveillance system or a waste heat recovery.


And of course you can not only equip it the way you like, but also design it the way you like. While our standard model is black on the outside and metallic on the inside, your personal LifeCube can come in any color you choose. Our color palette has eight different colors, which are freely available without extra charge. However, you can also pick other colors on request. In the end your LifeCube will fit perfectly to you, your studio, and your customers.



**inovatürk**

# ADVANTAGES OF COSMETIC-SOLUTION

Cosmetic-Solution manufactures all of its products in Germany. This way we can guarantee a consistent, excellent quality. Thanks to our twelve years of experience in the cosmetic sector, our products are technically mature and perfected. Yet, our service does not end with manufacturing or selling our products – even after we sealed the deal, we stay in contact with our clients and advise them in any way possible.



In this context, we do not only manufacture the LifeCubes ourselves, but also deliver and install them ourselves. Furthermore, we take care of any possible repairs, carry out the maintenances, and are only one phone call away if any questions arise.

At Cosmetic-Solution you do not have to buy a pig in a poke. You can visit us any time you want, examine our showroom, and even try out our LifeCube. Just make an appointment with us and we will be happy to show you around, discuss financing options, or give you further information about our products.

One part of our comprehensive service is a starter package, which includes several bathrobes, towels, slippers, headbands, and cotton gloves. This service is free of charge for all our LifeCube buyers. It enables you to start into business right away and welcome first customers immediately after the instalment of your personal cryo chamber.

To be well prepared for the first rush, we also offer individual training courses. They will take place at your location directly after the assemblage of the machine. During the instruction, we will explain the correct handling of the LifeCube as well as the interactive app to you.

*Get off to a flying start with cryotherapy  
and benefit from our electrical cryo chamber "LifeCube"*



**inovaturk**